

# BEHGAM Product Catalogue 2022

(Liners and Sleeves)



## Iceross Comfort Locking Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-5406XX	Transtibial	1 to 2	3mm
I-5006XX	Transtibial	1 to 2	6mm

Liner Size Standard: 16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40, 45.





LOCKING



STABILISING  
MATRIX



FABRIC  
COVER



## Iceross Comfort Cushion Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-CL53XX	Transtibial	1 to 2	3mm
I-CL56XX	Transtibial	1 to 2	6mm

Liner Size Standard: 16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40.



Size 40 is only available in 3mm.



CUSHION





FABRIC  
COVER



## Iceross Dermo Locking Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-4013XX	Transtibial	1 to 2	3mm
I-4016XX	Transtibial	1 to 2	6mm



Liner Size Standard:16,18,20,22,23.5,25,26.5,28,30,32,34,36,40,45.



## Iceross Dermo Cushion Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-CL43XX	Transtibial	1 to 3	3mm
I-CL46XX	Transtibial	1 to 3	6mm

Liner Size Standard:16,18,20,22,23.5,25,26.5,28,30,32,34,36,40.  
Size 40 is only available in 3mm.






## Iceross Seal-In X5 Transfemoral Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-TF673XX	Transfemoral	1 to 4	Standard 3mm
I-TF678XX	Transfemoral	1 to 4	Conical 3mm

Liner Size Standard:25,26.5,28,30,32,34,36,38,40,45,50,55.

Liner Size Conical:25,26.5,28,30,32,34,36,38.




## Iceross Seal-In Transfemoral Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-7532XX	Transfemoral	1 to 4	Standard 3mm
I-7632XX	Transfemoral	1 to 4	Conical 3mm

Liner Size Standard:25,26.5,28,30,32,34,36,38,40,45,50,55.

Liner Size Conical:25,26.5,28,30,32,34,36,38.





## Iceross Transfemoral Locking Liner (Ossur)

Article No	Amputation Level	Impact Level	Liner Profile
I-7032XX	Transfemoral	1 to 3	Standard 2mm
I-7132XX	Transfemoral	1 to 3	Conical 3mm

Liner Size Standard: 25, 26.5, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55.

Liner Size Conical: 25, 26.5, 28, 30, 32, 34, 36, 38.



## General Purpose Liner GPDT/GPFR (ALPS)

Article No	Suspension Method	Impact Level	Liner Profile
GPDTXX	Locking	1 to 2	3,6 mm
GPFRXX	Cushion	1 to 2	3,6 mm

Liner Size: 16,20,24,26,28,32,38,44.



LOCKING




CUSHION



GRIP GEL



BEIGE FABRIC  
WITH BOTH  
HORIZONTAL & VERTICAL  
STRETCH



## Extreme Liner AKDTHD/AKFRHD (ALPS)

Article No	Suspension Method	Impact Level	Liner Profile
AKDT(HD)XX	Locking	1 to 3	3,6 mm
AKFR(HD)XX	Cushion	1 to 3	3,6 mm

Liner Size: 16,20,24,26,28,32,38,44.



LOCKING



CUSHION




GRIP GEL



HIGH  
DENSITY GEL



BROWN FABRIC  
WITH REDUCED  
VERTICAL STRETCH





## Smart Seal Liner SAFR (ALPS)

Article No	Suspension Method	Impact Level	Liner Profile
SAFRXX	Cushion	1 to 4	3mm

Liner Size:26,28,32,38.



CUSHION



HYPOBARIC  
SUSPENSION



HIGH  
DENSITY GEL



IN HOUSE  
KNITTED FABRIC  
WITH HIGH  
QUALITY YARNS



## Easy Liner ELDT/ELFR (ALPS)

Article No	Suspension Method	Impact Level	Liner Profile
ELDTXX	Locking	1 to 4	3,6 mm
ELFRXX	Cushion	1 to 4	3,6 mm

Liner Size:16,20,24,26,28,32,38,44.



LOCKING



CUSHION




EASY GEL



SKIN FRIENDLY



BEIGE FABRIC  
WITH BOTH  
HORIZONTAL & VERTICAL  
STRETCH



## Superior Performance Liner SPDTHD/SPFRHD (ALPS)

Article No	Suspension Method	Impact Level	Liner Profile
SPDT(HD)XX	Locking	2 to 4	3,6 mm
SPFR(HD)XX	Cushion	2 to 4	3,6 mm



Liner Size: 16,20,24,26,28,32,38,44.



LOCKING



CUSHION



HIGH DENSITY GEL



GRIP GEL



BLACK FABRIC WITH HIGHLY REDUCED VERTICAL STRETCH



Suspension System



## LITE Liner Gel Excelsior (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile
4W7000XX	Locking	2 to 4	6mm
4W7500XX	Cushion	2 to 4	6mm

Liner Size:IS,M,L,IL,XL,XLP,XXL.



## LITE Liner Gel Prime (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile
4W8100XX	Locking	1 to 3	6mm
4W8600XX	Cushion	1 to 3	6mm

Liner Size:IS,M,L,IL,XL,XLP,XXL.



## LITE Liner Silicone SkinTex 2 TF Locking (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile	Amputation Level
4W2021XX	Locking	2 to 4	3mm	Transfemoral


Liner Size:26.5,28,30,32,34,36,38,40,42,45,50,55.



## LITE Liner Silicone Clear (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile	Amputation Level
4W1021XX	Locking	2 to 4	3mm	Transtibial
4W1121XX	Cushion	2 to 4	3mm	Transtibial

Liner Size:16,18,20,21,22,23.5,25,26.5,28,30,32,34,36,38,40,42,45,50,55.



## LITE Liner Silicone SkinTex 2 Soft (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile	Amputation Level
4W4022XX	Locking	1 to 3	3mm	Transtibial
4W4122XX	Cushion	1 to 3	3mm	Transtibial



Liner Size:18,20,21,22,23.5,25,26.5,28,30,32,34,36,38,40,42,45,50,55.

## LITE Liner Silicone SkinTex 6 Soft (Wagner)

Article No	Suspension Method	Impact Level	Liner Profile	Amputation Level
4W4062XX	Locking	1 to 3	6mm	Transtibial
4W4162XX	Cushion	1 to 3	6mm	Transtibial



Liner Size:18,20,21,22,23.5,25,26.5,28,30,32,34,36,38,40,42,45.



## Silicone Pads (Wagner)

Article No	Distal Gel Thickness	Diameter
4W15-08	12.5mm	80mm
4W15-10	11mm	100mm
4W15-12	10.7mm	120mm
4W15-14	12.8mm	140mm



## Distal Cup Silicone (Wagner)

Article No	Distal Thickness	Distal Circumference	for Silicone Socket Size
4W2000-06	5mm	6cm	14
4W2000-07	5mm	7cm	16
4W2000-09	5mm	9cm	18
4W2000-11	5mm	11cm	20/21
4W2000-12	5mm	12cm	22
4W2000-14	6mm	14cm	23.5
4W2000-16	6mm	16cm	25
4W2000-17	6mm	17cm	26.5
4W2000-19	6.5mm	19cm	28
4W2000-20	7mm	20cm	30
4W2000-22	7mm	22cm	32
4W2000-24	7mm	24cm	34



## Suspension Sleeve Gel with fabric cover (Wagner)

Article No	Total Length	Lower Circumference	Upper Circumference
4W80-XS	350mm	10-17cm	16-24cm
4W80-S	350mm	15-22cm	23-34cm
4W80-M	350mm	21-30cm	29-43cm
4W80-MP	350mm	24-35cm	33-47cm
4W80-L	350mm	28-42cm	36-51cm
4W80-LP	350mm	31-46cm	39-55cm
4W80-XL	350mm	34-50cm	42-60cm
4W80-XXL	350mm	36-55cm	44-70cm



Gel thickness: 3 mm.

Impact Level: 1 to 3.

## Antioxidant Cream (ALPS)

Article No	Amount
BP440A	50ml





# شرکت بهگام توان آرای سپاهان

تولید کننده و وارد کننده تجهیزات ارتز و پروتز

☎ 031-31314050

🌐 [www.behgamco.com](http://www.behgamco.com)

📷 [behgam\\_co](https://www.instagram.com/behgam_co)