


# BEHGAM Product Catalogue 2022

## (Lower Limb Prosthesis)



## Solid Ankle Foot with Toes (SACH Foot)

| Article No | Side | Size    | Thread | Weight Limit |
|------------|------|---------|--------|--------------|
| BP307B     | L/R  | 20~21cm | M10    | 90kg         |
| BP307B     | L/R  | 22~25cm | M10    | 95kg         |
| BP307B     | L/R  | 26~30cm | M10    | 100kg        |



## Solid Ankle Foot with Toes. Negroid (SACH Foot)

| Article No | Side | Size    | Thread | Weight Limit |
|------------|------|---------|--------|--------------|
| BP307B-1   | L/R  | 20~21cm | M10    | 90kg         |
| BP307B-1   | L/R  | 22~25cm | M10    | 95kg         |
| BP307B-1   | L/R  | 26~30cm | M10    | 100kg        |



## LP Solid Ankle Foot (SACH Foot)

| Article No | Size    | Thread | Weight Limit |
|------------|---------|--------|--------------|
| BP307LPB   | 21~27cm | M10    | 100kg        |



## Single Axis Foot with Two-hole

| Article No | Side | Size    | Weight Limit |
|------------|------|---------|--------------|
| BP309B     | L/R  | 21~25cm | 80kg         |
| BP309B     | L/R  | 26~30cm | 90kg         |



### SACH Foot for Women (Wagner)

| Article No | Side | Size    | Thread | Weight Limit |
|------------|------|---------|--------|--------------|
| BP307L25W  | L/R  | 21-26cm | M8     | 100kg        |



### SACH Foot for Men (Wagner)

| Article No | Side | Size    | Thread | Weight Limit |
|------------|------|---------|--------|--------------|
| BP307W     | L/R  | 22-29cm | M10    | 100kg        |



### Single Axis Foot for Men (Wagner)

| Article No | Side | Size    | Weight Limit |
|------------|------|---------|--------------|
| BP309W     | L/R  | 22-29cm | 100kg        |



## Pirogoff Foot Wagner for Men

| Article No | Side | Size    | Weight Limit |
|------------|------|---------|--------------|
| BP308WW    | L/R  | 22-29cm | 100kg        |



## Single Axis Foot Wagner With Blocks

| Article No | Side | Size    | Weight Limit |
|------------|------|---------|--------------|
| BP309WW    | L/R  | 21-28cm | 100kg        |



## Pirogoff Foot Without Blocks

| Article No | Side | Size    | Weight Limit |
|------------|------|---------|--------------|
| BP308B     | L/R  | 23-27cm | 90kg         |



## Footshell

| Article No | Side | Size    |
|------------|------|---------|
| BP346P     | L/R  | 22-30cm |



## Footshell

| Article No | Side | Size    |
|------------|------|---------|
| BP347P     | L/R  | 22-27cm |



## Footshell

| Article No | Side | Size    |
|------------|------|---------|
| BP348P     | L/R  | 22-30cm |



## Footshell

| Article No | Side | Size    |
|------------|------|---------|
| BP349P     | L/R  | 22-27cm |



## VARI-FLEX Foot

| Article No | Size    | Weight | Build Height | Heel Height | Weight Limit |
|------------|---------|--------|--------------|-------------|--------------|
| BP438CVO   | 22-30cm | 700g   | 170mm        | 10mm & 19mm | 166kg        |



## LP VARI-FLEX Foot

| Article No | Size    | Weight | Build Height | Heel Height | Weight Limit |
|------------|---------|--------|--------------|-------------|--------------|
| BP438CLPVO | 22-30cm | 680g   | 68mm         | 166kg       |              |



## FLEX-SYMES Foot

| Article No | Size    | Weight | Build Height | Heel Height | Weight Limit |
|------------|---------|--------|--------------|-------------|--------------|
| BP438CSO   | 22-30cm | 602g   | 60mm         | 10mm        | 166kg        |



## SACH Foot Adapter M10

| Article No | Material  | Weight | Weight Limit |
|------------|-----------|--------|--------------|
| BP310T10B  | Titanium  | 70g    | 125kg        |
| BP310S10B  | Steel     | 125g   | 125kg        |
| BP310A10B  | Aluminium | 80g    | 100kg        |



## Single Axis Foot Adapter

| Article No | Material | Size    | Weight | Weight Limit |
|------------|----------|---------|--------|--------------|
| BP311TSB   | Titanium | 22-25cm | 200g   | 100kg        |
| BP311TLB   | Titanium | 26-30cm | 210g   | 100kg        |
| BP311SSB   | Steel    | 22-25cm | 325g   | 100kg        |
| BP311SLB   | Steel    | 26-30cm | 340g   | 100kg        |



## Single Axis Foot Adapter Bumper

| Article No | Material |
|------------|----------|
| BP434FB    | Rubber   |



## Spectra Sock

| Article No | Size    |
|------------|---------|
| BP435S     | 22-25cm |
| BP435S     | 26-30cm |



## Connection plate for SACH Foot

| Article No | Material |
|------------|----------|
| BP312PSB   | Plastic  |



## Connection plate for Single Axis Foot

| Article No | Size    | Material |
|------------|---------|----------|
| BP312PSAB  | 22-29cm | Plastic  |



## Shaped Ankle Blocks

| Article No | Feet-Size      | Material | Bolt Assembly |
|------------|----------------|----------|---------------|
| BP437WB    | Men<br>22-29cm | Wood     | M10           |



## Bolt and Screw for SACH Foot

| Article No | Feet-Size      | Material              | Bolt Assembly |
|------------|----------------|-----------------------|---------------|
| BP329B     | Men<br>22-29cm | Steel<br>&<br>Plastic | M10           |





### Modular Single Axis Knee Joint with Lock (3R17)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP720SLB   | Steel    | 695g   | 150kg        |



### Modular Single Axis Knee Joint with Lock (3R33)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP320TLB   | Titanium | 530g   | 150kg        |



### Modular Polycentric Knee Joint (3R20)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP317SB    | Steel    | 625g   | 100kg        |



### Modular Polycentric Knee Joint (3R36)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP415TB    | Titanium | 445g   | 100kg        |



## Modular Disarticulation Knee Joint (3R21)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP315SB    | Steel    | 1010g  | 125kg        |



## Modular Knee Joint (3R30)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP413TB    | Titanium | 592g   | 125kg        |



## Modular Disarticulation Knee Joint with Lock (3R23)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP412SLB   | Steel    | 880g   | 125kg        |



## Modular Disarticulation Knee Joint with Lock (3R32)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP414TLB   | Titanium | 644g   | 125kg        |



## Modular Lightweight Single Axis Knee Joint with Lock (3R40)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP820ALB   | Aluminum | 286g   | 100kg        |



## Pyramid Adapter Wagner

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP409W     | Steel    | 95g    | 150kg        |
| BP410W     | Titanium | 68g    | 150kg        |



## Modular Hip Joint (7E7)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP313EB    | Titanium | 685g   | 100kg        |



## Modular Monocentric Knee Joint with Friction Brake (3R15)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP319SB    | Steel    | 490g   | 100kg        |



## Modular Monocentric Knee Joint with Friction Brake (3R49)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP416TB    | Titanium | 360g   | 100kg        |



### Modular Safety Knee Joint Wagner (3R15)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP319SW    | Steel    | 650g   | 100kg        |



### Modular Safety Knee Joint Wagner (3R49)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP416TW    | Titanium | 501g   | 100kg        |



### Modular Knee Joint Wagner (3R17)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP720SLW   | Steel    | 560g   | 125kg        |



### Modular Knee Joint Wagner (3R33)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP320TLW   | Titanium | 419g   | 125kg        |



### Modular Knee Joint Wagner (3R23)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP412SLW   | Steel    | 939g   | 125kg        |



### Modular Knee Joint Wagner (3R32)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP412TLW   | Titanium | 644g   | 125kg        |



### Modular Lightweight Knee Joint Wagner (3R40)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP820ALW   | Aluminum | 279g   | 85kg         |



### Modular Knee Joint Wagner (3R20)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP317SLW   | Steel    | 701g   | 100kg        |



### Modular Knee Joint Wagner (3R36)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP415TW    | Titanium | 462g   | 100kg        |



### Modular Knee Joint Wagner (3R21)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP315SW    | Steel    | 887g   | 125kg        |



### Modular Knee Joint Wagner (3R30)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP413TW    | Titanium | 592g   | 125kg        |



## Rotation Adapter

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP411SB    | Steel    | 298g   | 100kg        |
| BP411TB    | Titanium | -      | 100kg        |



## Rotation Adapter (Wagner)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP411SW    | Steel    | 295g   | 100kg        |
| BP411TW    | Titanium | 190g   | 100kg        |



## Screw Rotation Adapter (Wagner)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP411TTW   | Titanium | 236g   | 125kg        |
| BP412TTW   | Titanium | 263g   | 125kg        |



**3R21 Anchor**



**7E7 Anchor**



**Extension Assist (3R15)**



**Extension Assist (3R21)**



**Extension Assist Cable (3R15)**



**System Manual Lock (3R23)**





## Tube Clamp Adapter 10 Degree

| Article No | Material | Diameter | Weight | Weight Limit |
|------------|----------|----------|--------|--------------|
| BP431SAB   | Steel    | 30mm     | 133g   | 100kg        |
| BP431TAB   | Titanium | 30mm     | 85g    | 100kg        |



## Tube Clamp Adapters

| Article No | Material | Diameter | Weight | Weight Limit |
|------------|----------|----------|--------|--------------|
| BP402SB    | Steel    | 30mm     | 132g   | 100kg        |
| BP402AB    | Aluminum | 30mm     | 75g    | 100kg        |
| BP402TB    | Titanium | 30mm     | 86g    | 100kg        |
| BP402S34B  | Steel    | 34mm     | 140g   | 150kg        |
| BP402T34B  | Titanium | 34mm     | 85g    | 150kg        |



## Tube Clamp Adapters, Sliding

| Article No | Material | Diameter | Weight | Weight Limit |
|------------|----------|----------|--------|--------------|
| BP323SSB   | Steel    | 30mm     | 281g   | 100kg        |
| BP323ASB   | Aluminum | 30mm     | 156g   | 75kg         |
| BP323TSB   | Titanium | 34mm     | 185g   | 100kg        |



## Tube Adapters

| Article No | Lenght | Material | Diameter | Weight | Weight Limit |
|------------|--------|----------|----------|--------|--------------|
| BP333SL34B | 40cm   | Steel    | 34mm     | 370g   | 150kg        |
| BP333TL34B | 40cm   | Titanium | 34mm     | 330g   | 150kg        |



## Tube Adapters

| Article No | Lenght | Material | Diameter | Weight | Weight Limit |
|------------|--------|----------|----------|--------|--------------|
| BP333S40B  | 40cm   | Steel    | 30mm     | 295g   | 100kg        |
| BP333S25SB | 20cm   | Steel    | 30mm     | 199g   | 100kg        |
| BP333A40B  | 40cm   | Aluminum | 30mm     | 275g   | 100kg        |
| BP333A25B  | 20cm   | Aluminum | 30mm     | 179g   | 100kg        |
| BP333T40B  | 40cm   | Titanium | 30mm     | 267g   | 100kg        |
| BP333T25B  | 20cm   | Titanium | 30mm     | 170g   | 100kg        |
| BP333C40B  | 40cm   | Carbon   | 30mm     | 117g   | 100kg        |
| BP333C25B  | 20cm   | Carbon   | 30mm     | 58g    | 100kg        |



## Extension Adapter With Pyramid Receiver

| Article No | Lenght | Material | Diameter | Weight | Weight Limit |
|------------|--------|----------|----------|--------|--------------|
| BP333S60B  | 60mm   | Steel    | 30mm     | 126g   | 100kg        |
| BP333S90B  | 90mm   | Steel    | 30mm     | 160g   | 100kg        |
| BP333T60B  | 60mm   | Titanium | 30mm     | -      | 100kg        |
| BP333T90B  | 90mm   | Titanium | 30mm     | -      | 100kg        |



## Tube Adapters (Wagner)

| Article No | Material | Lenght | Weight | Weight Limit |
|------------|----------|--------|--------|--------------|
| BP333T40W  | Titanium | 40cm   | 271g   | 125kg        |
| BP333T20W  | Titanium | 20cm   | 171g   | 125kg        |
| BP333S40W  | Steel    | 40cm   | 323g   | 125kg        |
| BP333S20W  | Steel    | 20cm   | 221g   | 125kg        |



## Tube Clamp Adapters (Wagner)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP402TW    | Titanium | 73g    | 125kg        |
| BP402SW    | Steel    | 116g   | 125kg        |
| BP402AW    | Aluminum | 72g    | 100kg        |



## Double Adapter Slidable (60mm)

| Article No | Material | Sliding | Weight | Weight Limit |
|------------|----------|---------|--------|--------------|
| BP323TS60B | Titanium | ±11mm   | 208g   | 100kg        |
| BP323SS60B | Steel    | ±11mm   | 330g   | 100kg        |



## Double Pyramid & 4 Screw Adapters

| Article No | Lenght | Material | Weight | Weight Limit |
|------------|--------|----------|--------|--------------|
| BP432SFM   | 75mm   | Steel    | 185g   | 125kg        |
| BP432SFM   | 60mm   | Steel    | 160g   | 125kg        |
| BP432SFM   | 45mm   | Steel    | 135g   | 125kg        |
| BP432SFM   | 32mm   | Steel    | 110g   | 125kg        |
| BP432TFM   | 75mm   | Titanium | 115g   | 100kg        |
| BP432TFM   | 60mm   | Titanium | 98g    | 100kg        |
| BP432TFM   | 45mm   | Titanium | 88g    | 100kg        |
| BP432TFM   | 32mm   | Titanium | 72g    | 100kg        |



## Double Pyramid Adapters

| Article No | Lenght | Material | Weight | Weight Limit |
|------------|--------|----------|--------|--------------|
| BP433SMM   | 32mm   | Steel    | 115g   | 125kg        |
| BP433TMM   | 32mm   | Titanium | 73g    | 125kg        |



## Double Receiver Adapters

| Article No | Lenght | Material | Weight | Weight Limit |
|------------|--------|----------|--------|--------------|
| BP326S75B  | 75mm   | Steel    | 205g   | 125kg        |
| BP326S60B  | 60mm   | Steel    | 185g   | 125kg        |
| BP326S45B  | 45mm   | Steel    | 162g   | 125kg        |
| BP326S32B  | 32mm   | Steel    | 116g   | 125kg        |
| BP326A75B  | 75mm   | Aluminum | 106g   | 100kg        |
| BP326A60B  | 60mm   | Aluminum | 92g    | 100kg        |
| BP326A45B  | 45mm   | Aluminum | 87g    | 100kg        |
| BP326A32B  | 32mm   | Aluminum | 69g    | 100kg        |
| BP326T75B  | 75mm   | Titanium | 125g   | 125kg        |
| BP326T60B  | 60mm   | Titanium | 110g   | 125kg        |
| BP326T45B  | 45mm   | Titanium | 95g    | 125kg        |
| BP326T32B  | 32mm   | Titanium | 85g    | 125kg        |



## Pyramid Adapter with hole

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP322SSB   | Steel    | 95g    | 100kg        |
| BP322TSB   | Titanium | 60g    | 100kg        |



## Rotatable Pyramid Adapter with 4 holes

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP322STB   | Steel    | 160g   | 125kg        |



## Socket Adapter with 4 holes

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP321SSB   | Steel    | 120g   | 100kg        |



## Rotatable Socket Adapter with 4 holes

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP321STB   | Steel    | 160g   | 125kg        |



## Base of Pyramid Adapter

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP417STB   | Steel    | 95g    | 125kg        |



## 3 Prong Pyramid Adapter

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP325SSB   | Steel    | 125g   | 100kg        |



## 3 Prong Rotatable Socket Adapter with Pyramid

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP325STB   | Steel    | 175g   | 125kg        |
| BP325TTB   | Titanium | 140g   | 125kg        |



## 3 Prong Rotatable Socket Adapter with Receiver

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP419STB   | Steel    | 175g   | 125kg        |
| BP419TTB   | Titanium | 145g   | 125kg        |



## Rotatable Socket Adapter with 4 ears

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP328STB   | Steel    | 150g   | 100kg        |



## Socket Adapter with 4 ears

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP328SSB   | Steel    | 105g   | 100kg        |



## Lamination Anchor with 4 prongs

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP324SH    | Steel    | 105g   | 125kg        |





## Pyramid Adapter

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP418SSB   | Steel    | 70g    | 100kg        |
| BP418TSB   | Titanium | 40g    | 100kg        |



## Receiver Adapter

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP420SSB   | Steel    | 70g    | 100kg        |
| BP420TSB   | Titanium | 50g    | 100kg        |



## Shuttle Lock with Pyramid Adapter

| Article No | Material | Distal Connection | Weight Limit |
|------------|----------|-------------------|--------------|
| BP336MS    | Aluminum | Pyramid Adapter   | 100kg        |



## Shuttle Lock for Silicone Liner

| Article No | Material | Distal Connection | Weight Limit |
|------------|----------|-------------------|--------------|
| BP337FS    | Aluminum | Adjustment Screw  | 100kg        |



## Shuttle Lock Below Knee

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP336MB    | Aluminum | 301g   | 100kg        |



## Shuttle Lock Above Knee

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP337FB    | Aluminum | 304g   | 100kg        |



## Plastic Shuttle Lock

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP338PB    | Plastic  | 156g   | 100kg        |



### Shuttle Lock Pyramide (Wagner)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP336MW    | Aluminum | 191g   | 100kg        |



### Shuttle Lock Receiver (Wagner)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP337FW    | Aluminum | 201g   | 100kg        |



### Shuttle Lock (Wagner)

| Article No | Material | Mounting Height | Weight |
|------------|----------|-----------------|--------|
| BP338PW    | Plastic  | 37mm            | 99g    |



### Shuttle Lock (ALPS)

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP338PA    | Plastic  | 126g   | 100kg        |



## Screw Valve

| Article No | Material | Weight Limit |
|------------|----------|--------------|
| BP425PB    | Aluminum | 100kg        |



## Flat Rubber Valve

| Article No | Outside Ø of Seat Ring | for Hole Ø |
|------------|------------------------|------------|
| BP427LB    | 40mm                   | 24mm       |
| BP427LB    | without seat ring      | 24mm       |



## Listra Flat Combi Valve Wagner (Automatic)

| Article No | Outside Ø of Seat Ring | for Hole Ø |
|------------|------------------------|------------|
| BP428VW    | 40mm                   | 24mm       |
| BP428VW    | without seat ring      | 24mm       |



## Flat Silicone Valve

| Article No | Material | Weight Limit |
|------------|----------|--------------|
| BP426SB    | Silicone | 100kg        |



## Icelock Expulsion Valve Ossur

| Article No | Material | Weight Limit |
|------------|----------|--------------|
| BP428FO    | Plastic  | 100kg        |



## Socket Attachment Block

| Article No | Material | Weight | Weight Limit |
|------------|----------|--------|--------------|
| BP404WB    | PU-ABS   | 275g   | 100kg        |



## Connection Tube with Seat Ring (Wagner)

| Article No | Outside Ø of Seat Ring | for Hole Ø |
|------------|------------------------|------------|
| BF423AKW   | 40mm                   | 24mm       |



## Connection Tube with Seat Ring

| Article No | Material | Weight Limit |
|------------|----------|--------------|
| BP423PB    | Plastic  | 100kg        |



## Perlon Transfemoral Cosmetic Stockings - AK (Wagner)

| Article No | Leg Length | Cuff Circumference | Foot Length |
|------------|------------|--------------------|-------------|
| BM139AKW   | 55cm       | 41cm               | 20cm        |
| BM139AKW   | 65cm       | 41cm               | 20cm        |
| BM139AKW   | 75cm       | 43cm               | 20.5cm      |
| BM139AKW   | 80cm       | 45cm               | 21cm        |
| BM139AKW   | 85cm       | 45cm               | 22cm        |



## Perlon Transtibial Cosmetic Stockings - BK (Wagner)

| Article No | Leg Length | Cuff Circumference | Foot Length |
|------------|------------|--------------------|-------------|
| BM139BKW   | 36cm       | 25cm               | 20cm        |
| BM139BKW   | 41cm       | 25cm               | 21cm        |
| BM139BKW   | 46cm       | 27cm               | 22cm        |



## Perlon Cosmetic Stockings

| Article No | Amputation Level |
|------------|------------------|
| BM139AKB   | Transfemoral     |
| BM139BKB   | Transtibial      |



## Advanced Transfemoral Foam Cover

| Article No | Lenght | Material          | Calf Circumference | Knee Flexion | Side |
|------------|--------|-------------------|--------------------|--------------|------|
| BM105AKAB  | 80cm   | PUR flexible foam | 36cm,40cm,44cm     | 35°          | L/R  |



## Transfemoral Foam Cover

| Article No | Lenght | Material          | Calf Circumference | Knee Flexion | Side |
|------------|--------|-------------------|--------------------|--------------|------|
| BM105AKB   | 80cm   | PUR flexible foam | 36cm,40cm,44cm     | 30°          | L/R  |



## Transfemoral Foam Cover (Reticular)

| Article No | Size | Material | Dimension |
|------------|------|----------|-----------|
| BM105AKHB  | L    | EVA foam | 18×85cm   |
| BM105AKHB  | M    | EVA foam | 17.5×83cm |
| BM105AKHB  | S    | EVA foam | 17×80cm   |



## Transtibial Foam Cover (Plastazote)

| Article No | Lenght | Material | for Tube Diameter |
|------------|--------|----------|-------------------|
| BM105PBKB  | 48cm   | PE foam  | 30cm              |



## Transtibial Foam Cover

| Article No | Lenght | Material          | for Tube Diameter |
|------------|--------|-------------------|-------------------|
| BM105BKB   | 55cm   | PUR flexible foam | 30cm,34cm         |







# شرکت بهگام توان آرای سپاهان

تولید کننده و وارد کننده تجهیزات ارتز و پروتز

☎ 031-31314050

🌐 [www.behgamco.com](http://www.behgamco.com)

📷 [behgam\\_co](https://www.instagram.com/behgam_co)