


BEHGAM Product Catalogue 2022

(Upper Limb Prosthesis)



Upper Limb



8E38=8 Sensor Hand Speed (Ottobock)

Article No	Side	Size	Inner Hand
8E38=8-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E38=8-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E38=8-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E38=8-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E38=8-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E38=8-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Reference No	8E38=8	8E38=8	8E38=8
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7.2 V	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	15-300 mm/sec	15-300 mm/sec	15-300 mm/sec
Proportional speed	0-100 N	0-100 N	0-100 N
Weight with system inner hand	462 g	462 g	462 g
for	Women Adolescents	Men	Men



8K22 Single-Cable System Hand (Ottobock)

Article No	Side	Threaded stud	Size	Inner hand
8K22=L6 3/4	Left	M12×1.5	6 3/4	8X18=L6 3/4
8K22=L7 1/4	Left	M12×1.5	7 1/4	8X14=L7 1/4
8K22=L7 3/4	Left	M12×1.5	7 3/4	8X14=L7 3/4
8K22=L8	Left	M12×1.5	8	8X14=L8
8K22=R6 3/4	Right	M12×1.5	6 3/4	8X14=R6 3/4
8K22=R7 1/4	Right	M12×1.5	7 1/4	8X14=R7 1/4
8K22=R7 3/4	Right	M12×1.5	7 3/4	8X14=R7 3/4
8K22=R8	Right	M12×1.5	8	8X14=R8

Reference Number	8K22	8K22	8K22	8K22
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 215 g	ca. 300 g	ca. 330 g	ca. 340 g



8S11 Prosthetic Glove for Adolescents and Men (Ottobock)

Article No	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S11=190X76L	Left	7 1/4	8X18=L7 1/4 8X24=L7 1/4	300 mm	230 mm	Adolescents Men
8S11=190X76R	Right	7 1/4	8X18=R7 1/4 8X24=R7 1/4	300 mm	230 mm	Adolescents Men
8S11=192X78L	Left	7 1/4	8X18=L7 1/4 8X24=L7 1/4	440 mm	245 mm	Adolescents Men
8S11=192X78R	Right	7 1/4	8X24=R7 1/4	440 mm	245 mm	Adolescents Men
8S11=210X78L	Left	7 3/4	8X18=L7 3/4 8X24=R7 3/4	320 mm	250 mm	Men
8S11=210X78R	Right	7 3/4	8X18=R7 3/4 8X24=R7 3/4	320 mm	250 mm	Men
8S4=220X80L	Left	8	8X14=L8	225 mm	210 mm	Men
8S4=220X80R	Right	8	8X14=R8	225 mm	210 mm	Men
8S11=225X80L	Left	8 1/4	8X18=L8 1/4 8X24=L8 1/4	350 mm	260 mm	Men
8S11=225X80R	Right	8 1/4	8X24=R8 1/4	350 mm	260 mm	Men



8S6=170x65 Prosthetic Glove for Children (Ottobock)

Article No	For hand size	Side	Sleeve length
8S6=170x65L	6 3/4	Left	300 mm
8S6=170x65R	6 3/4	Right	300 mm



8S5=195x78 Prosthetic Glove for Women (Ottobock)

Article No	For hand size	Side	Sleeve length
8S5=195x78L	7 1/4	Left	340 mm
8S5=195x78R	7 1/4	Right	340 mm



8S4=210x78 Prosthetic Glove for Men (Ottobock)

Article No	For hand size	Side	Sleeve length
8S4=210x78L	7 3/4	Left	220 mm
8S4=210x78R	7 3/4	Right	220 mm



Female Soft Voluntary Opening (SVO) Hand (Fillauer) - Left

Article No	Size	Wrist to Thumb	Cable Exit	Thread ½-20	Thread M12	Glove Number
BP103CF	7 ¼ in.	5 ¼ in.	Dorsal	61064	61997	20701
BP103CF	7 ¼ in.	5 ¼ in.	Palmar	61484	61999	20701

Female Soft Voluntary Opening (SVO) Hand (Fillauer) - Right

Article No	Size	Wrist to Thumb	Cable Exit	Thread ½-20	Thread M12	Glove Number
BP103CF	7 ¼ in.	5 ¼ in.	Dorsal	61063	61996	20702
BP103CF	7 ¼ in.	5 ¼ in.	Palmar	61483	61998	20702



Male Soft Voluntary Opening (SVO) Hand (Fillauer) - Left

Article No	Size	Wrist to Thumb	Cable Exit	Thread ½-20	Thread M12	Glove Number
BP104CF	7 ¾ in.	5 ⅜ in.	Standard	61066	61993	30703
BP104CF	8 ¼ in.	5 ⅜ in.	Standard	61068	62006	30707
BP104CF	7 ¾ in.	5 ⅜ in.	Palmar	61486	61991	30703
BP104CF	8 ¼ in.	5 ⅜ in.	Palmar	61488	62008	30707


Male Soft Voluntary Opening (SVO) Hand (Fillauer) - Right

Article No	Size	Wrist to Thumb	Cable Exit	Thread ½-20	Thread M12	Glove Number
BP104CF	7 ¾ in.	5 ⅜ in.	Standard	61065	61994	30704
BP104CF	8 ¼ in.	5 ⅜ in.	Standard	61067	62005	30708
BP104CF	7 ¾ in.	5 ⅜ in.	Palmar	61485	61992	30704
BP104CF	8 ¼ in.	5 ⅜ in.	Palmar	61487	62007	30708



Gloves for Child Passive Hands (Fillauer)

Article No	Side	Approximate Age for Use	A	B	C	D	E
10900L 10900R	Left Right	0-18months	100mm 100mm	100mm 100mm	40mm 40mm	133mm 133mm	175mm 175mm
10901 10902	Left Right	1-3years	115mm 115mm	115mm 115mm	45mm 45mm	145mm 145mm	210mm 210mm
10903 10904	Left Right	3-6years	140mm 140mm	115mm 115mm	50mm 50mm	145mm 150mm	230mm 230mm
10905 10906	Left Right	5-9years	150mm 150mm	130mm 130mm	55mm 55mm	160mm 165mm	255mm 255mm



10900L
10900R



10901
10902



10903
10904



10905
10906

Upper Limb


نحوه انتخاب دست پسیو:

دست پسیو = ۰۰۰ + کد کالا + ۵

مثال: کد کالا ۱۰۹۰۱ است، شماره دست پسیو، ۱۰۰۰۰۰۹۰۱۰۵ می‌شود.

Gloves for Female Passive Hands (Fillauer)

Article No	Side	Remarks	A	B	C	D	E
22907	Left	Medium nails	154mm	122mm	70mm	200mm	410mm
22908	Right		151mm	123mm	65mm	200mm	400mm
22909	Left	Medium nails	167mm	128mm	72mm	214mm	440mm
22910	Right		165mm	133mm	70mm	215mm	440mm
22911	Left	Long nails	170mm	148mm	75mm	205mm	490mm
22912	Right		170mm	148mm	75mm	217mm	490mm
22703	Left	Long nails	174mm	140mm	78mm	222mm	470mm
22704	Right		172mm	138mm	78mm	222mm	470mm
22917	Left	Long nails Prominent veins	180mm	157mm	82mm	245mm	430mm
22918	Right		184mm	157mm	82mm	245mm	430mm
20909	Left	Unisex Short nails	162mm	134mm	75mm	195mm	320mm
20910	Right		164mm	142mm	73mm	204mm	320mm
20911	Left	Unisex Short nails	170mm	148mm	75mm	210mm	380mm
20912	Right		170mm	148mm	75mm	235mm	380mm
20913	Left	Long Fingers Long nails	172mm	146mm	90mm	215mm	370mm
20914	Right		170mm	142mm	90mm	223mm	370mm
20915	Left	Short nails	175mm	150mm	84mm	220mm	360mm
20916	Right		178mm	150mm	82mm	199mm	380mm
20917	Left	Prominent veins	190mm	162mm	90mm	230mm	385mm
20918	Right		198mm	158mm	85mm	200mm	385mm



22907 22909 22911 22703 22917 20909 20911 20913 20915 20917
22908 22910 22912 22704 22918 20910 20912 20914 20916 20918

نحوه انتخاب دست پسیو:


دست پسیو = ۰۰۰ + کد کال + ۵

مثال: کد کال ۲۲۹۰۷ است، شماره دست پسیو، ۵۲۲۹۰۷۰۰۰ می‌شود.

Upper Limb

Gloves for Male Passive Hands (Fillauer)

Article No	Side	Remarks	A	B	C	D	E
30909	Left		171mm	143mm	75mm	206mm	330mm
30910	Right		168mm	146mm	75mm	215mm	330mm
30911	Left		185mm	156mm	83mm	215mm	320mm
30912	Right		185mm	156mm	83mm	210mm	320mm
30913	Left	Prominent	197mm	167mm	90mm	230mm	370mm
30914	Right	Veins	200mm	161mm	85mm	225mm <td 360mm	
30915	Left		190mm	146mm	82mm	235mm	370mm
30916	Right		200mm	168mm	90mm	240mm	450mm
30917	Left	Long Fingers	214mm	178mm	95mm	224mm	400mm
30918	Right		225mm	186mm	100mm	295mm	400mm
30919	Left		217mm	176mm	92mm	275mm	360mm
30920	Right		212mm	174mm	90mm	240mm	340mm
30921	Left		222mm	186mm	95mm	282mm	360mm
30922	Right		225mm	177mm	95mm	285mm	370mm



30909
30910

30911
30912

30913
30914

30915
30916

30917
30918

30919
30920

30921
30922


نحوه انتخاب دست پسیو:

دست پسیو = ۰۰۰ + کد کالا + ۵

مثال: کد کالا ۳۰۹۰۹ است، شماره دست پسیو، ۵۳۰۹۰۹۰۰۰ می‌شود.


Female Myo/Mechanical Glove

Article No	Size	A	B	C
BP105CB	Large	175mm	190mm	75mm
BP105CB	Medium	170mm	180mm	70mm
BP105CB	Small	145mm	170mm	65mm



Male Myo/Mechanical Glove

Article No	Size	A	B	C
BP105MCB	Large	180mm	195mm	80mm
BP105MCB	Medium	175mm	190mm	75mm
BP105MCB	Small	175mm	180mm	70mm




Silicone Finger Set

Article No	Gender	Size	Length
BP105FB	Female	Large	80mm
BP105FB		Medium	75mm
BP105FB		Small	68-70mm
BP106FB	Male	Large	85mm
BP106FB		Medium	80-78mm
BP106FB		Small	70mm




Female Silicone Cosmetic Glove

Article No	Size	A	B	C
BP105B	Large	160mm	195mm	80mm
BP105B	Medium	155mm	180mm	75mm
BP105B	Small	140mm	165mm	68-70mm




Male Silicone Cosmetic Glove

Article No	Size	A	B	C
BP105MB	Large	180mm	220mm	85mm
BP105MB	Medium	165mm	200mm	80-78mm
BP105MB	Small	155mm	185mm	70mm




Child Silicone Cosmetic Glove

Article No	Years	A	B	C
BP105BB	2-3	115mm	115mm	45mm
BP105BB	3-6	115mm	140mm	50mm
BP105BB	5-9	130mm	150mm	55mm




Female Cosmetic Glove with Zip

Article No	Size	A	B	C
BP110F	Large	160mm	195mm	80mm
BP110F	Medium	155mm	180mm	75mm
BP110F	Small	140mm	165mm	68-70mm




Male Cosmetic Glove with Zip

Article No	Size	A	B	C
BP110M	Large	180mm	220mm	85mm
BP110M	Medium	165mm	200mm	80-78mm
BP110M	Small	155mm	185mm	70mm




Female Cosmetic Glove with Inner Filler and Bolt

Article No	Size	A	B	C
BP111F	Large	160mm	195mm	80mm
BP111F	Medium	155mm	180mm	75mm
BP111F	Small	140mm	165mm	68-70mm




Male Cosmetic Glove with Inner Filler and Bolt

Article No	Size	A	B	C
BP111M	Large	180mm	220mm	85mm
BP111M	Medium	165mm	200mm	80-78mm
BP111M	Small	155mm	185mm	70mm



Mechanical Hand (ALNA)

Article No	Threaded stud	Size	Weight
BP109	M12×1.5	7 3/4	185 g



Upper Limb

E-200 Elbow (Hosmer)

Article No	Color	Cable Exit	Weight	Locking Positions	Saddle	Gear Sector
BP113	Tan	Outside	350g	11	Standard	Standard



WE Friction Wrist

Article No	Lamination Ring	Material	Diameter	Weight	Thread
BP116L45B	No	Aluminium	44mm	57g	½-20
BP116L45B	Yes	Aluminium	44mm	65g	½-20
BP116L50B	No	Aluminium	51mm	74g	½-20
BP116L50B	Yes	Aluminium	51mm	78g	½-20



Wooden Wrist Fitting

Article No	Diameter	Thread
BP115B	60mm	M12×1.5mm




Elbow Set-Up

Article No	Upper Arm Lamination Ring	Wrist Inner
BP112LB	67mm	45mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R810	80mm	110mm	225mm	128mm	240mm	191mm	271mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L810	80mm	110mm	220mm	127mm	240mm	184mm	265mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R811	71mm	105mm	210mm	111mm	240mm	176mm	242mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L811	71mm	105mm	214mm	112mm	240mm	184mm	245mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R812	80mm	110mm	225mm	123mm	240mm	188mm	269mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L812	80mm	110mm	225mm	122mm	240mm	188mm	265mm




Gloves for Male and Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R813	75mm	106mm	198mm	120mm	240mm	165mm	238mm




Gloves for Male and Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L813	75mm	106mm	198mm	122mm	240mm	169mm	240mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R814	84mm	115mm	225mm	125mm	210mm	185mm	260mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L814	84mm	115mm	225mm	122mm	210mm	194mm	266mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R815	73mm	105mm	212mm	120mm	240mm	184mm	256mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L815	73mm	115mm	210mm	117mm	240mm	180mm	257mm




Gloves for Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R816	72mm	100mm	183mm	111mm	210mm	154mm	223mm




Gloves for Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L816	72mm	100mm	182mm	110mm	210mm	155mm	221mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R817	88mm	110mm	215mm	128mm	240mm	179mm	252mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L817	88mm	110mm	210mm	128mm	240mm	179mm	253mm




Gloves for Children Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R818	70mm	90mm	173mm	105mm	190mm	142mm	205mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L818	70mm	90mm	165mm	103mm	190mm	142mm	198mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R819	80mm	110mm	217mm	127mm	240mm	175mm	245mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L819	80mm	110mm	212mm	125mm	240mm	175mm	247mm




Gloves for Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R820	74mm	105mm	187mm	119mm	210mm	164mm	238mm




Gloves for Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L820	74mm	105mm	189mm	115mm	210mm	168mm	238mm




Gloves for Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R821	80mm	110mm	194mm	119mm	220mm	167mm	241mm




Gloves for Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L821	80mm	110mm	192mm	116mm	220mm	162mm	235mm




Gloves for Male and Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R822	70mm	105mm	200mm	116mm	220mm	166mm	237mm




Gloves for Male and Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L822	70mm	105mm	204mm	104mm	240mm	167mm	235mm




Gloves for Male and Female Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R823	92mm	105mm	218mm	121mm	220mm	170mm	237mm




Gloves for Male and Female Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L823	80mm	110mm	207mm	127mm	220mm	172mm	235mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R824	80mm	105mm	220mm	118mm	240mm	183mm	253mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L824	80mm	110mm	210mm	123mm	240mm	178mm	247mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R825	75mm	105mm	215mm	118mm	240mm	178mm	242mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L825	75mm	105mm	213mm	115mm	240mm	180mm	239mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R826	70mm	100mm	192mm	115mm	240mm	172mm	237mm




Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L826	70mm	100mm	192mm	113mm	240mm	173mm	232mm




Gloves for Male Total Hands (Right)

Article No	C1	C2	F	D	B	A0	A18
R828	85mm	115mm	222mm	126mm	240mm	183mm	252mm



Gloves for Male Total Hands (Left)

Article No	C1	C2	F	D	B	A0	A18
L828	85mm	115mm	222mm	126mm	240mm	184mm	254mm





شرکت بهگام توان آرای سپاهان

تولید کننده و وارد کننده تجهیزات ارتز و پروتز

☎ 031-31314050

🌐 www.behgamco.com

📷 [behgam_co](https://www.instagram.com/behgam_co)